

Rules in the Cell lab.

- The door must always be closed
- The window may not be open

- Always clean with 70 % EtOH before using the LAF bench.
- Always clean with both soap (deconex) and 70 % EtOH after using the LAF bench.
- Always clean under the plates in the LAF bench.

- Remember to turn off LAF bench, microscope, centrifuge (leave lid for centrifuge open), and water bath after use.

- Remember to cover up microscopes.

- You may not change and save in the programs already made at the centrifuge, but you may create new.

- At the end of the day, turn on the UV light in all LAF.

- In general, tidy up after use.

- Remember to write name and date on everything you are using.

- If you use the last, or nearly the last of the flasks with 70 % ETOH or Deconex, you must fill it up.
- If you use the last of the 5 liter 70 % ETOH, get a new. You can find it in the chemistry room A401-118.

- Remember to fill up in the cupboards and drawers. You will find more material in the "BRIK" cupboards.

- The sterilized glassware you use in the cell lab. You can find them in the sterilization oven in room A401-166. After use, rinse the glassware in d- water and place the glassware in the dishwasher tub in Lab A401-104.

- If you wish to have items sterilized, put them on top of the oven, and we will see to, that they are sterilized overnight, and will be ready the next day.

- You can use the cell-counter, but you must buy your own counting chamber slide.

- You must follow the rules for working in a GMO I Lab. You can read the specifications on the door.